



RED LEAF



MANUAL SUCTION PUMP

SP-01-Ultra


Technical features

Hand held, hand powered, portable suction unit with overflow safe guard and container for aspirated contents which can easily be disassembled to be washed and reassembled and re-used.
Includes lumen reducing adaptor and carry case.

Suction power adjustable
50% For Child
100% For Adult



Fr 12 Catheter connection
Exter Diameter: Φ 4mm



Fr 14 Catheter connection
Exter Diameter: Φ 4.67mm



Child Catheter connection
Exter Diameter: Φ 9mm



Adult Catheter connection
Exter Diameter: Φ 17mm

Model	Unfolded Size (L x W x H)	Folded Size (L x W x H)	Self-weight	Load Bearing	Packing		
					Quantity	Size (L x W x H)	G.W
SP-01	-	-	-	-	25	50x40x45cm	14kg