

ACB-264SL (ZVC-264SL) (Vaccine Carrier – Long Range)

Model Reference	ACB-264SL (ZVC-264SL)
Type	Long Range Carrier
Hot Zone	YES
Cold Chain Compliant	YES

Specifications

Gross Storage Capacity (L)	11.6
Vaccine Storage Capacity (L)	6
External Surface Material	LLDPE*
Internal Lining Material	LLDPE*
Insulation Material	CFC free PU
Insulation Thickness	100mm

Lid Type & Fixings	Hinged
Number of Ice Packs Required	26
Number of Ice Packs Supplied	26
Ice Pack Type	E5/IP.2
Volume of Ice Pack (L)	0.40

Performance (E04/CB01)

Cold Life without opening	132 Hours @ 43°C
---------------------------	------------------

Dimensions

Boxed Dimensions (mm)	(L) 650mm x (W) 530mm x (H) 460mm
External Dimensions (mm)	(L) 630mm x (W) 510mm x (H) 450mm
Internal Dimensions (mm)	(L) 400mm x (W) 280mm x (H) 206mm
Vaccine Storage Dimensions (mm)	(L) 330mm x (W) 210mm x (H) 100mm
Weight (Kg's) - Fully Loaded	26.6
Weight (Kg's) - empty (with empty icepacks)	14.21
Weight Boxed (Kg's) - empty (with empty icepacks)	16



Supplied By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

🌐 albert@zeromedical.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted