

## MED690R (Single Swing Glass Door Medical Fridge with Low - E Glass)

|                     |     |
|---------------------|-----|
| ISO 9001            | YES |
| IEC 60335 Compliant | YES |

### Specifications

|   |  |
|---|--|
| Manufacturers Gross Volume (L)            | Refrigerator 525L  |
| Adjustable Shelves                        | 4 (Stainless Steel)  |
| Inner Body Material                       | White Coated Mild Steel (Chromadek)  |
| Exterior Body Material                    | White Coated Mild Steel (Chromadek)  |
| Fridge Temperature Display                | YES  |
| Temperature (Fridge)                      | +2°C To +8°C At +35°C Ambient  |
| Accurate Temperature Control              | YES (0.1°C Accuracy)   |
| Interior Light - LED                      | YES (Auto Off After 10 Min Of Operating Door)  |
| Visible Digital Temperature Control       | YES  |
| Internal Compartment Fan                  | YES  |
| Door Locking System                       | YES  |
| Glass Doors                               | 1 – Double Glazed With UV Protection (Low-E)   |
| Automatic Defrost                         | YES  |
| High / Low Temperature Alarm +2°C to +8°C | YES (Audible & Visual & Will Sound If Temperature Out Of Set Parameters For More Than 5 Minutes) |
| Door Open Alarm                           | YES (Audible Alarm Will Sound After 60 Seconds When The Door Has Been Left Open)                 |
| Battery Alarm For Power Outage            | YES (Audible & Visual With >48Hrs Capacity)  |
| Insulation                                | Ecomate Polyurethane (Zero GWP & Zero ODP)   |
| Door Open Alarm                           | YES (Visual After 60 Seconds)  |
| CFC & HCFC Free                           | YES  |
| Mains Cable Length With Plug              | 3.1M   |
| Adjustable Feet                           | Optional   |
| Swivel Castors                            | YES - 4 (2 Lockable)   |
| Voltage Stabiliser                        | YES (Combined Surge Protect & Voltage Stabiliser)  |
| RTMD Ready                                | YES  |
| Interior Light Manual Switch              | YES (External)   |



### Performance Information

|                                 |              |
|---------------------------------|--------------|
| Hold Over Time During Power Cut | 4 To 8 Hours |
|---------------------------------|--------------|

### Dimensions

|                          |                                    |
|--------------------------|------------------------------------|
| External Without Castors | (H) 2020mm x (W) 690mm x (D) 635mm |
| External With Castors    | (H) 2120mm x (W) 690mm x (D) 635mm |
| Weight (Kg's) – Boxed    | 120                                |



Supplied By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

🌐 albert@zeromedical.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted