



## BATTERY ASPIRATOR

### 478-000 - Port-a-Vac

#### Technical features

Fully portable battery operated and made from durable plastic shell with 12V5AH re-chargeable battery inside, and can be used with vehicle power supply DC 12V or charge by the city electric AC220V. (All connections included) Battery life indicator as well as suction gauge are incorporated into the unit.

Suction bottle is clearly marked with level indicator and suction gauge.

Has over flow cut off valve.

Two year warranty.



#### Included accessories:

- Extreme pressure bearing value: 650mmHg
- Pressure bearing indicator scope: 0- 760mmHg
- Vacuum pressure: -50~-750mmHg
- Air extracting rate:25L/min
- Noise: ≤60dB(A)
- Liquid storage bottle capacity: 800ml
- Suction tubing: 10mm x 1.5m
- Power supply: AC220V±10%, 50Hz or DC12V
- Self-weight: 5kg
- Size:35x28x16cm

