

Examination light, LED, Mobile

KL01L1LED Suitable for ENT, urology, gynecology examination and outpatient surgery half auxiliary lighting.

- F-features**
- ◆ An imported LED cold light source is adopted as operation lighting. As a genuine cold light source.
 - ◆ The LED shadowless lamp is pure-DC power supplied, ensuring no flickering or dimming occurs.
 - ◆ The service life of LED shadowless lamp is long (50,000 h).
 - ◆ LED has high luminous efficiency and it is impact-resistant and strong to resist crumbling, without mercury contamination. The light contains no infra-red or ultraviolet radiation pollution.
 - ◆ The unique design of the spring arm. Allows for easy, single-handed manoeuvring.
 - ◆ A weighted base allows for manoevrability without compromising on balance.
 - ◆ Effective lighting over the bed with a minimum length of 1 metre.

Technical Specifications

The light may be separately switched on via the column or on the luminaire.

Colour rendering index: > 90

Working distance of 50 cm = 15 cm

Working distance of 100 cm = 25 cm

Operates from 230V, 50 Hz mains power

Mobile stand: lockable durable castors + hook for mains cable

Smooth surface allows for easy cleaning

Colour temperature 4000k

Lamp life:>50000 Continuous hours

Base diameter 40cm

60000 lux



Sterilisable positioning handle