

## MED85F (Vaccine Freezer Electric - Compression)



Version	MK1
Application	Fridge
Door Types	Solid Door (SD) Glass Door (GD)
ISO 9001	YES
IEC 60335 Compliant	YES

### Specifications

Manufacturers Nett Volume (L)	83L
Manufacturers Gross Volume (L)	93L
Adjustable Shelves	3
Number Of Vaccine Boxes (100mm x 100mm x 100mm)	= 60
Number Of Vaccine Boxes (100mm x 100mm x 50mm)	= 10
Compartment Light	YES (LED)
Inner Body Material	Aluminium
Exterior Body Material	Stainless Steel & White Coated Mild Steel (Chromadek)
Freezer Temperature Display	YES
Visible Digital Temperature Control	YES
Internal Compartment Fan	YES
Door Locking System	YES
Glass Door Specification	Double Glazed With UV Protection (Low-E)
Door Open Alarm (Audible)	YES (Audible Alarm Will Sound After 60 Seconds When The Door Has Been Left Open)
Door Open Alarm (Visual)	YES (Visual After 60 Seconds)
CFC & HCFC Free	YES
Mains Cable Length With Plug	2.3M
Adjustable Feet	YES
Interior Light Manual Switch	YES (External)
Swivel Castors	NO (Optional, Add 100mm To Height)

### Performance Information

Solid Door Unit	-15°C To -25°C
Glass Door Unit	-15°C To -20°C

### Dimensions

Internal	(H) 530mm x (W) 420mm x (D) 420mm
Boxed Dimensions	(H) 985mm x (W) 640mm x (D) 605mm
Unboxed Dimensions	(H) 975mm x (W) 620mm x (D) 585mm
Weight (Kg's) - Boxed	61

### Container Quantities

20'	54
40'	114
40' HC	114



Manufactured By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

🌐 albert@zeromedical.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted