


10 LITRE OXYGEN CONCENTRATOR

Model 8F-10



- Automatic reminder filter cleaning: 100 hours
- Flow Range: 1~10L/min
- Oxygen Concentration: 95.5%~87%
- Maximum Limited Pressure: 100kPa
- Sound Level: 55dB(A)
- Input Power: 700VA
- Minimum Operating Time: 30 minutes



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
8F-10	1	23.5	41x37.5x66.8