



SIRIO BABY 200

Portable ICU ventilator for children and newborns

Simple ventilation advanced Technology

The project of the new **Sirio Baby 200** started in 1982 and now, this friendly user equipment has been refined achieving the goal of combining an exceptional flexibility with a very high level technology.

This ventilator has been conceived especially for the emergency when a ventilatory support is needed during the transport. Solid and light at the same time and equipped with a long duration back-up battery which grants the maximum flexibility.

Versatility

The new Sirio Baby 200 is able to offer a wide range of therapeutic treatments as well as to perform a pressure and volume control ventilation guaranteed either with an invasive or with a non-invasive modality, able to meet all the respiratory needs of newborns and prematures (weight < 3kgs).

Simplicity

The intuitive and versatile interface carefully studied for doctors operating with neonatal emergency situations makes this ventilator the ideal equipment to be used along with the transport incubator and for the transport in ambulance, mobile clinics, helicopter and airplanes.

Performances

The new Sirio Baby 200 grants excellent performances taking advantage of Siare's high standards technology.

Equipped with different ventilation modalities: CPAP – IMV – IPPV – MANUAL – PEEP – Pressure trigger – FiO2 adjustable from 21% to 100%. Inspiratory and Expiratory time adjustable. Integrated spirometry.

Long duration back-up battery and a wide range of dedicated accessories grant an easier mobility of the device.

MRI Version

The new Sirio Baby 200 is manufactured also in the amagnetic version (item code 09637/RMN) for the use in the Magnetic Resonance unit (MRI) from 1,5 to 3,0 Tesla. The amagnetic version is the ideal choice for the transport of newborns from the resuscitation room to the magnetic resonance unit with a continue ventilation support granting to the newborns a higher comfort and security.

SIRIO BABY 200





The new Sirio Baby 200 can be used in Magnetic Resonance (MRI) from 1,5 to 3,0 Tesla.



The new Sirio Baby 200 can be used also in helicopter, mobile clinics, ambulances and airplanes.

Acoustic and visual alarms

Simple and intuitive settings

Built-in spirometer



Adjustable trigger

Long duration battery

Patient circuit with expiratory valve and proximal flow sensor



Transport version with professional bag for ventilator and 2 It cylinder with inox plate for being hooked/unhooked in the ambulance.



Transport version with professional inox support for ventilator and 2 lt cylinder with support for the relevant installation in the ambulance.

by Siare M

DATA FOR ORDERING






Item Code 09637 SIRIO BABY 200 Neonatal Portable ICU Ventilator



Item Code 09637/RMN SIRIO BABY 200 Neonatal Portable ICU Ventilator for MRI

OPTIONAL ACCESSORIES



For additional accessories please refer to the price list in force

Distributed by:







