



## N5A Nasal Mask

A balance of breathing and horizon

Forehead frame –free design enables you to enjoy more vision.



**AutoFit Cushion**

Smart groove design on the cushion helps to automatically and effectively adjust on the nose bridge.

**Compact Design**

N5A is composed of 3 parts only. Easy to use and clean.

**Liberty of Vision**

Based on the forehead frame-free design, users can read before they go to sleep.

**Easy Release Clips**

Easy release clips provide the user with more ease of use.

**Product Codes**

**N5A Nasal Mask**

- 10020193 (L)
- 10020194 (M)
- 10020195 (S)



**Frame Assembly**

60380278



**Cushion Assembly**

- 60380242 (L)
- 60380243 (M)
- 60380244 (S)



**Headgear Assembly**

61030052



**BMC Medical Co., Ltd.**