DATA FOR ORDERING



Code 09436-SIRIO S2/T



Code A26.049048 Stand for ventilator with 5 wheels

Code G00208001

with housing for

O₂ cylinder

complete with

for accessories.

ventilator and 2 lt.

safety band and bag

Stainless steel support



Code G00203001 Plate for rapid connection of tubolar support 10G (code G00208001)

on wall.



Code 102641/SLU

subject to change without

Silicone adult/paediatric patient circuit with expiratory valve. A70.049001





Code G00203000

Plate for rapid connection of ventilator on ambulance flat



Code 10006/2

Pin Index O2 cylinder 2 It capacity, empty

Code 109006/CM

Pin Index O2 pressure reducer



Code G90003010

Padded and refracting transport bag for ventilator and 2lt O2 cylinder complete with ambulance connection.







Via G. Pastore, 18 - Località Crespellano 40053 Valsamoggia (BO) Italy Phone: +39 051 969802 (r.a.) - Fax: +39 051 969366 E N G I N E E R I N G INTERNATIONAL GROUP s.r.l. General e-mail: mail@siare.it - Internet: www.siare.it











EMERGENCY

SIRIO S2/T

Emergency lung ventilator



Easy ventilation, Advanced technology

The Sirio S2/T, a project started in 1979 and a user-friendly equipment that combines excellent flexibility and advanced technology.

It has been expressly conceived for patient's breathing assistance in emergency healthcare during transport. Compact and light, it is equipped with a long lasting internal battery granting the maximum flexibility of use.

Versatility:

Sirio S2/T offers a wide range of therapeutic options; it is able to ensure a controlled volume ventilation limited by pressure in invasive and non-invasive modality, meeting the breathing requirements of adult, paediatric and neonatal patients (weight > 3 kg).

Simplicity:

The user-friendly interface, cleverly designed to facilitate emergency health operators, makes of this ventilator the ideal choice for use on ambulance, helicopter and airplanes.

Performances:

Sirio S2/T guarantees excellent performances thanks to the technological advantages granted by Siare.

It is equipped with several ventilation modes: CPAP - PSV (Apnoea Back-up) - IPPV+AST - MANUAL - PEEP - Pressure Trigger - FiO2 adjustable at 50% or 100%.

Long lasting internal battery and a series of optional accessories designed for a good portability of the device.

Sirio S2/T Professionalism in Emergency





Version for transport with professional bag for ventilator and 2 lt cylinder housing, complete with support for rapid connection/disconnection for installation in ambulance.



Version for transport with professional stainless-steel tubular support for ventilator and 2/3 lt cylinder housing, with possibility to attach a quick connection plate for installation in ambulance.

Certified according to EN 1789 Directive.



The Sirio S2/T is the ideal solution for patient transport.



The Sirio S2/T is homologated for use in helicopter.