

yuwell

OXYGEN CONCENTRATOR

8F-5A



Global sales leader



Medical grade



High oxygen concentration



5L oxygen flow



Failure detection



Ultra-silence



Castors

- Flow Range: 0.5~5L /min
- Oxygen Concentration: 95.5%~87%
- Maximum Limited Pressure: 70kPa
- Sound Level: 5L/min Sound Pressure ≤ 54dB(A)
- Input Power: 350VA
- Minimum Operating Time: 15minutes



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
8F-5A	1	19/16.5	47.5×33×58.5

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Manufacturing Address: No.1 Baisheng Road, Development Zone, Danyang,

Jiangsu Province, P.R.China

Tel: 86-0511-8690 0859 86-0511-8690 0860

E-mail: yuwellsales@yuyue.com.cn; sales@yuyue.com.cn

www.yuwell.com

f | t | in | i | Yuwell Global