



Fix and Free

iVolve F1A Full Face Mask

BMC-FM1A

Superior performance delivers maximum comfort. With lighter design and effective treatment, FM1A will offer you and your family a good night sleep.



Ultra comfort

Dual-wall cushion will significantly fit for your face with minimal pressure points. The slide-in design provides you with high level seal and effective treatment.

Friendly usability

Thanks to its 5 parts compact design, you can reassemble it instantly. Less parts, more convenience.

Visual Freedom

With the unique design, your vision is totally released and you will have the chance to enjoy the world as well as therapy.

Product Codes

Thin Pad & Cushion

S (61030016)

M (61030017)

L (61030018)



Headgear Frame

(60350272)



Soft Sleeves

(60380173)



BMC Medical Co., Ltd.

Add: 5/F Main Building, No.19 Gucheng Street West, Shijingshan, Beijing 100043, China

Tel: +86-10-51663880

Fax: +86-10-51663880 ext.810

en.bmc-medical.com