

ZVFF240AC VACCINE STORAGE REFRIGERATOR AND ICE PACK FREEZER – ELECTRIC COMPRESSION

International Cold Chain Compliant	YES
ISO 9001 : 2015	YES
IEC 60335 Compliant	YES

Specifications

Dual Cooling System	Refrigerator (Separate On/Off Switch) Freezer (Separate On/Off Switch)
Manufacturers Gross Volume (L)	Refrigerator 205L Freezer 40L
Shelves	5 (Stainless Steel)
Inner Body Material	Stainless Steel
Fridge Temperature Display	YES
Freezer Temperature Display	YES
Solid Door Finish	YES
Temperature (Fridge)	+2°C To +8°C At +35°C Ambient
Temperature (Freezer)	-10°C To -15°C At +35°C Ambient
Accurate Temperature Control	YES
Hold Over Facility	YES
Separate Freezer Door	YES
Interior Light When Door Open	YES
Stainless Steel Inner Liner	YES
Visible Digital Temperature Control	YES (0.1°C Accuracy)
Internal Compartment Fan	YES
Door Locking System	YES
Door	With Internal Illumination & Control Panel
Protection Against Freezing The Product	YES
Automatic Defrost	YES
Visual Alarm System For Sensor Failure	YES
High / Low Temperature Alarm +2°C To +8°C	YES (Audible & Visual & Will Sound If Temperature Out Of Set Parameters For More Than 5 Minutes)
Voltage Stabiliser	YES (Combined Surge Protect & Voltage Stabiliser)
Door Open Alarm	YES (Audible Alarm Will Sound After 60 Seconds When The Door Has Been Left Open)
Battery Alarm For Power Outage	YES (Audible & Visual With >48Hrs Capacity)
CFC & HCFC Free	YES
Mains Cable Length With Plug	3.1M
Adjustable Feet	Optional (Remove 100mm From Height)
Swivel Castors	YES - 4 (2 Lockable)
Power Consumption (kWh)	0.050 (Stable Running At +35°C)
RTMD Ready	YES

Performance Information

Hold Over Time During Power Cut	>48 Hours (Depending On Ambient)
Ice Pack Freezing	8Kg/24Hrs

Dimensions

Boxed Dimensions (mm)	(H) 1930mm x (W) 675mm x (D) 660mm
Unboxed Dimensions (mm)	(H) 1900mm x (W) 630mm x (D) 620mm
Weight (Kg's) - boxed	115kg

HOLD OVER UNIT



Manufactured By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

✉ albert@zeromedical.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted