



# NET CYLINDERS

## NEW ENERGY TECHNOLOGY CO., LTD.

Portable medical oxygen respirators design and manufacture by NET, are produced through advanced technology. Comparing with conventional respirators, we use aluminium composite cylinders instead of steel ones, with weight reduction of 40%-70%. More preferable for people with breathing problems. Full specifications meet our customers' demands. NET supplies two types of O<sub>2</sub> cylinders for all your choice, all aluminium cylinders and aluminium composite cylinders.

### Technical Data - All Aluminium Cylinders

- Complies with International Safety Standards: EN 1975:1999+A1:2003 and EN ISO 7866:2012+AC:2014
- Aluminium Cylinder without built-in regulator
- Pin Index Yoke Fitting
- Bullnose Fitting
- Can be filled from a high pressure oxygen booster pump

