

# CHAMPCARE INTERNATIONAL LIMITED

29/F, Nanyang Plaza, 57 Hung To Road, Kwun Tong, Kowloon, Hong Kong  
Tel: (852) 3106 4048

## Ultra DGD25A Pneumatic (Air driven) Gas Booster Pump / System



---

**Portable 2-Stage Pneumatic (Air Driven) High Pressure Oxygen Refilling Booster system, including 200 litre Compressor with Purge Valve to release excessive high pressures.  
Oxygen is compressed up to 200 BAR.**

### **FEATURES:**

- Portable
- Pin Index Yoke to fit portable cylinders
- Couples to Bullnose outlet
- High Pressure Outlet (steel) hosing with Bullnose fitting to fit 50L Oxygen Cylinders (6m)
- High Pressure Outlet (steel) flexible hosing with Bullnose fitting (3m)
- Low Pressure hosing as attachment to compressor and booster pump (6m)
- Inlet, outlet and pressure gauges.

Champcare International Limited

Address: 29/F, Nanyang Plaza, 57 Hung To Road, Kwun Tong, Kowloon, Hong Kong

Tel. : 852 3106 4048

Email : [patrick.tai@champcare-limited.com](mailto:patrick.tai@champcare-limited.com)