








OXYGEN CONCENTRATOR

8F-10 SINGLE FLOW

			
Medical grade	High oxygen concentration	10L oxygen flow	Cumulative time function
			
Failure detection	Castors	High efficiency particulate air Filter	

- Automatic reminder filter cleaning: 100 hours
- Flow range: 1-10L/min
- Oxygen concentration : 95.5%-87%
- Maximum limited pressure : 100kPa
- Sound level : 55dB(A)
- Input power : 700VA
- Minimum operating time : 30minutes



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	N.W (KG)	Size (cm)
8F-10 SINGLE FLOW	1	23.5	41×37.5×66.8

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Sales Lenter Address: Huanyuan tast Road No. 1, Xuzhuang Software Park, Nanjing

Jianasu Province, PR. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn

www.yuwell.com