



FOETAL DOPPLER

BF-500D+ & BF-500D+ (TFT) Ultrasound Foetal Doppler

Features

- Portable design, the probe holder with a card-type protection seat design, easy loading and unloading, with a gallus, easy to carry;
- Probe with highly efficient interface design, stable performance, continuous wave transducer 2 Crystal narrow beam, 2.0MHz ultrasound probe to detect from 10 weeks;
- LCD displays 3-digit FHR accurately, syncheart rate flicker signal, battery symbol, volume and foetal heart rate abnormal sound and visual alarm
- Shut down automatically in one minute without input and have foetal heart rate over range alarm and battery power alarm
- BF-500D+ with backlight for night use, the backlight can be set by the user to save energy and extend battery life
- BF-500D+ TFT can set the interface, count mode (four mode) alarm on/off and foetal heart rate range alarm settings, language switching etc.
- Built-in speaker, with a headphone jack, can transmit audio signals
- Multiple frequency probes for options to meet different needs;
- Audio response: 300Hz - 1KHz
- Output response: 9M Wmm²
- Power output of obstetric probes: 10sW/cm²
- Battery 9V Alkaline, estimated battery life: 8 hours from fully charged
- Self weight: 280g
- Dimensions: 15cm x 7.5cm x 4cm
- On / Off Button
- Volume controlled button
- Displays heart rate and heart beat synchronous indicator
- Low battery: Battery charging
- 220V 60Hz
- Supplied in a dust free bag with zip
- 3 Year Warranty



Shenzhen Bestman Instrument Co.,Ltd.

Add:4/F, Block 210, 2nd Industrial Area of Nanyou, Xiangnan Rd.,
Nanshan Dist., Shenzhen, P.R.China
Tel: 0086-755-26411423/26416184
Fax: 0086-755-26416184/26415407
Email: szbestman@szbestman.com
bestman@szbestman.com
szbestman@163.com
Website: <http://www.szbestman.com>
<http://www.bestmanhk.com>



ISO 13485:2003



ISO 9001:2008



www.szbestman.com