

yuwell

BATTERY PORTABLE SUCTION UNIT

Model 7E-G1



Oil-free
pump



Suction
liner



Vehicle
use



Over-flow
protection



Lithium
battery



Wall-hung
shelf

- DC power supply: AC100-240V DC 12V
- Input power: 150VA
- Max Vacuum: ± 85 kPa
- Air Flow Rate: 27L ± 4 /min
- Operating noise: ≤ 65 dB (A)
- Operating time: 30 minutes on, 30 minutes off
- Running time of battery operation: Approx. 60mins



Ecomed
Medical (Pty) Ltd

www.ecomed.co.za |
homecare@ecomed.co.za

25 Clew Street, Monument Ext. 1, Krugersdorp, 1739
+27 11 955 5710 | sales@ecomed.co.za