



KL-LED·MSTZ4



KL-LED·MSTZ4(AC/DC)

Mobile Theatre Light

Specifications

Mode	KL-LED.MSTZ4
Illumination at 1m(lux)	80000~160000
Light field diameter at 1m (cm)	10-28
D10(mm)	180±40
Diameter of spot D50(mm)	110±35
Color rendering index(Ra)	100≥Ra≥85

Mode	KL-LED.MSTZ4
Optional color temperature(K)	(3500-4500)±500
Number of LED bulbs	72
Average bulb life (hours)	50,000
Heat-to-light ratio(mW/m ² lx)	<6
Working distance(cm)	70-140
Max. power consumption(VA)	90