

Product Specifications

These click style flowmeters are designed for use in ambulances, hospitals, and portable oxygen systems. The click style flowmeter is easy to use because you simply turn the dial to obtain accurate flow setting. Unlike a traditional flowmeter, the click style flowmeter does not need to be in an up-right position in order to function properly. This flowmeter is light weight, rugged, and has large, easy to read numbers.

- Calibrated inlet pressure: 50 PSIG / 70 Deg F (21 Deg C)
- Operating temperatures: 20°F to 130°F (-28.8°C to 54.4°C)
- Inlet port: 1/8 FNPT to connect O2 heyer probe
- Outlet port: 9/16-18 UNF (CGA 1240)
- Storage temperature: 14°F to 125°F (-10°C to 51.6°C)
- Available in Oxygen
- Oxygen is available in 25, 15, 8 and 4 LPM,
- Available in single and double style flowmeter
- Smallest and lightest in the industry
- Smallest increment is 1/32
- Color coded labelling for specific gas
- Aluminum body color anodized, available in black
- Easy to read window
- ABS plastic knob
- Brass fixed orifice disc inside



**Ultra Single or Double
Single Dial-a-Flow meter**

**Ultra Single or Double
Double Dial-a-Flow meter**

