



KL-LED·STZ4 / STZ4

Ceiling Mounted Theatre Light

Specifications

Diameter of lighthouse(cm)	Combined type	Number of LED bulbs	72/72
Illumination at 1 m(lux)	80000-160000/80000-160000	Average bulb life (hours)	50000
Light field diameter at 1m(cm)	10-28	Heat-to-light ratio(mW/m ² lx)	<6
D10(mm)	180±40/180±40	Working distance (cm)	70-140
Diameter of 'spot DS0(mm)	110±35/110±35	Power supply	100-240V 50/60HZ
L1+L2 Depth of illumination(cm)	150	Input power(VA)	180
Color rendering index(Ra)	100≥Ra≥85	Ceiling height for installation(cm)	270-320
Optional color temperature(K)	6700≥Tc≥3000(K)		