

# BREATHWEAR MASK SERIES

## YN-02



Fit



Comfort



Convenience



Light & Soft

- Type: Nasal Mask
- Therapy Pressure Range :4-30 cmH2O
- Tube Connection :  $\phi$ 22mm
- Flow Resistance[At 50 L/min]:  $\leq$ 0.5cmH2O
- Flow Resistance[At 100 L/min]:  $\leq$ 1.5cmH2O



## PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YN-02	10	1.35/0.95	35×28×28