

Diagnostic Equipment

BP Meter Aneroid Multicuff Set

Specifications:

- Palm type blood pressure monitor
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Chrome casing 65mm diameter dial
- Black powder-coated handle
- Navy blue washable cuffs. Sizes: Infant, Child, Adult, Obese and Thigh
- Single Tube latex bladder in each cuff
- Standard latex bulb and valve



Ecomed
Medical (Pty) Ltd

t: +27 11 955 5710
f: +27 11 955 6316

www.ecomed.co.za
sales@ecomed.co.za