

## Automatic Blood Pressure Monitor

Model : AV151 **f**



- Real Fuzzy Technology
- Irregular Heartbeat Detector (IHB)
- 90 memories
- One-touch automatic operation
- String bag
- Latex-free patented cone cuff



BHS A/A graded approved



ESH approved



Real Fuzzy Technology



Irregular Heartbeat Detector (IHB)



90 memories



One-touch operation



Latex-free patented Cone cuff



**Clinically Validated**

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements and achieved the highest possible grade (A/A) according to the criteria of the British Hypertension Society. This monitor further complies with the standard required by the Association for the Advancement of Medical Instrumentation (AAMI).

### Real **f**uzzy Technology

Real Fuzzy Technology determines ideal cuff pressure based on one's systolic blood pressure and arm size. Users no longer need to pre-set the inflation level before measurement. The technology eliminates arm discomfort through unnecessary high inflation settings and also prevents inaccurate readings caused by erroneous cuff inflation levels.



Cone cuff fits arm circumference:  
24-36cm (9.4" ~14.2")  
Other size cuffs available upon request

- Storage and Transportation Condition  
-10°C ~ 60°C (14°F ~ 140°F)  
10% ~ 90% RH