

ZLF30AC VACCINE STORAGE REFRIGERATOR ELECTRIC – COMPRESSION WITH SURECHILL® TECHNOLOGY




ZERO
 A P P L I A N C E S

WHO PQS Approved

WHO PQS no:	E003/051
WHO Test Procedure	E003/RF03-VP.2
WHO Specifications	E003/RF03.2

Specifications

Vaccine Storage Capacity (L)	Refrigerator	27L
Frozen Ice Pack Storage Capacity (L)	Freezer	N/A
Manufacturers Gross Volume (L)	Refrigerator	30L
	Freezer	N/A
Hot Zone Appliance (+43°C)		YES
Freeze Protection		Grade A
CCEOP Extended Range Product		YES



Power Consumption – Day/Night Cycle

(kW/24Hrs) on 220 – 240V mains Electric	+32°C	0.620
	+43°C	0.940

Performance Information

Hold Over Time During Power Cut At +43°C on The Day / Night Cycle	> 110 Hours
Hold Over Time During Power Cut At +32°C on The Day / Night Cycle	> 150 Hours
Ice Pack Freezing	N/A

*** **Sure Chill**® Technology guarantees that vaccines cannot and will not freeze within rated temperatures (+10°C to +43°C)

Dimensions

Boxed Dimensions (mm)	(H) 1050mm x (W) 635mm x (D) 595mm
Unboxed Dimensions (mm)	(H) 1028mm x (W) 619mm x (D) 563mm
Weight (Kg's) - boxed	65kg
Shipping Volume (Crated)	0.372 m³
Weight Crated	100kg

Container Quantities

20'	54
40'	114
40' HC	114

Zero Appliances (Pty) Ltd • Johannesburg, South Africa

 010 207 1600 •
  010 207 1675 •
 Int  +27 10 207 1600 •
 Int  +27 10 207 1675

 info@zeroappliances.co.za •
  www.zeroappliances.co.za

E&OE Errors and Omissions Excepted