

ZVC-46 (Vaccine Carrier – long range)

Model Reference	ZVC-46
Type	Long Range Carrier
WHO PQS	YES – Hot Zone

Specifications

Vaccine Storage Capacity (L)	2.5
Weight - fully loaded	6.36kg
Weight - empty (with empty icepacks)	2.98kg
External Surface Material	HDPE (High Density PE)
Internal Lining Material	HIPS (High Impact PS)
Insulation Material	CFC free PU
Insulation Thickness	32/35mm

Lid Type & Fixings	Removable
Number of Ice Packs Required	4
Number of Ice Packs Supplied	4
Ice Pack Type	E5/IP.1

Performance

Cold Life without opening	43°C	50 Hours
---------------------------	------	----------

Dimensions

External Dimensions (mm)	(L) 1270mm x (W) 1270mm x (H) 320mm
Internal Dimensions (mm)	(L) 1220mm x (W) 1220mm x (H) 190mm
Vaccine Storage Dimensions (mm) with frozen ice pack	(L) 113.8mm x (W) 113.8mm x (H) 190mm

