

# Diagnostic Equipment

## BP Meter Aneroid Pocket Set

### Specifications:

- Palm type blood pressure monitor
- Deluxe plastic casing with 65mm diameter dial
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Black powder-coated handle
- Navy blue or red washable adult size nylon cuff
- 1 or 2-tube adult size cuff and bladder
- Standard latex bulb and valve
- Vinyl zipper case



**Ecomed**  
Medical (Pty) Ltd

t: +27 11 955 5710  
f: +27 11 955 6316

[www.ecomed.co.za](http://www.ecomed.co.za)  
[sales@ecomed.co.za](mailto:sales@ecomed.co.za)