Diagnostic Equipment

BP Meter Aneroid Square Face Wall Mounted with Swivel

Specifications:

- Square shaped 145 mm diameter dial
- Front and rear adjustable for easy reading
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- PVC coiled tubing





t: +27 11 955 5710 f: +27 11 955 6316 www.ecomed.co.za sales@ecomed.co.za