

# PORT-A-VAC

- Dimensions: Approx 34cm x 28cm x 16cm
- Weight: 4.5kg
- Variable Regulator to Adjust Vacuum Level
- Vacuum Range 150mmHg - 750mmHg (584mmHg  $\pm$ 76mmHg)
- Meets Ambulance Specifications for flow of 25Lpm for 60 minutes at full battery
- Protected Vacuum Gauge
- Reusable Collection Canister
- Bacteria Filter (2) supplied
- DC/Battery Power Options - Operates from Internal 12V Rechargeable Battery
- Exterior Charging Indicator Light
- Battery Charging Transformer Supplied
- 12V Vehicle Plug Supplied
- Battery automatically charges when plugged into 12V DC car lighter-outlet or 220V AC socket using our Transformer 478-251 (supplied)
- One Year Limited Warranty on Unit (Disposables Excluded)

