

# BF- series Ultrasonic Fetal Doppler



Note:BF-500+/500++/500B FDA approval

## Features:

- Pocket BF-500: It has no display;
- Pocket BF-500+:LCD can display FHR accurately;
- Pocket BF-500++:Blue backlight LCD displays FHR accurately, convenient to use at night;
- BF-500+/500++: Alarm function for abnormal FHR;
- High quality loudspeaker with clear sound and the audio output interface;
- Volume digital control;
- With headphone jack, Ni-MH rechargeable battery;
- With a gullus, easy to carry;
- Normative collocation: 2.5 MHz probe, easy to replace.



BF-500/500+/500++



BF-500A/500B

## Features:

- Pocket BF-500A:LCD can display FHR, synchronous FHR signal, battery capacity indicator, volume and abnormal alarm signal etc;
- Built-in rechargeable battery. Modem inside can dial up and transmit the fetal heart beat signal to the tele-monitor center through a telephone line;
- Pocket BF-500B: the same with BF-500A excluding the dial function;
- Automatic power off (one minute) and alarm function for abnormal heart rate;
- Collocated with gullus and portable bag, easy to carry;
- Normative collocation: 2.0 MHz probe, 5.0 MHz probe to detect umbilical blood flow (for option).