

MED1360R (Glass Door Medical Cooler With Low E Glass)

Model Reference

MED1360R

Features

White sides / black trim	Standard
Stainless Steel Sides/ grey trim	Optional
LED Lights	Standard
Electronic Controller (Thermostat With Display)	Standard
Over/Under Voltage Compressor Protection	Standard
Lock	Standard
Front Ventilation	Standard

Specifications

Gross Capacity	Ltr	1080
Nett Capacity	Ltr	900
Shelves	8 (extra shelves optional)	
Refrigerant	R134a	
Temp Range	2°C to 8°C	
Max Power Consumption (watt)	780	
Voltage / Hertz	230/50	
Compressor Size (hp/cc)	1/2/18	
Weight Empty (kg)	208	

Dimensions

External	(H) 2020mm x (W) 1360mm x (D) 635mm
Internal (Nett)	(H) 1380mm x (W) 1270mm x (D) 550mm



medical refrigeration
MEDCO



Medco Medical Consulting (Pty) Ltd • Johannesburg, South Africa

☎ 011 849 1758 • 📠 011 849 1758 •

🌐 info@medicalconsulting.co.za • 🌐 www.medicalconsulting.co.za

E&OE Errors and Omissions Excepted