

# Regulator - Pin Index High Flow Heyer



Dial Click Stop Constant Flow with  
pressure gauge to register cylinder pressure  
Oxygen delivery through Click-Stop  
Hand Dial Marked 0-25l/min  
Oxygen Cylinder Pressure Gauge  
Pressure Dial protected by a rubber ring  
Outer Body manufactured from Aluminium  
Inner Body Manufactured from Brass  
Minimum Input Operating Pressure 200 bars  
Output Pressure at Heyer: 3 - 6 Bar  
Male/Female Connection Included

2 Year Warranty & 3x Bodoc Seals included