



**KL01L.II LED**



● KL01L.II LED (AC/DC)

● KL01L.II LED



## KL01L.II LED

*KL01L.IILED Suitable for ENT, urology, gynecology examination and outpatient surgery half auxiliary lighting.*

**F-  
eatures**

- ◆ An imported LED cold light source is adopted as operation lighting. As a genuine cold light source.
- ◆ The LED shadowless lamp is pure-DC power supplied, without any strobe or harmonic.
- ◆ The service life of LED shadowless lamp is long (50,000 h).
- ◆ LED has high luminous efficiency and it is impact-resistant and strong to resist crumbling, without mercury contamination. The light it sends contains no infra-red or ultraviolet radiation pollution.
- ◆ The unique design of the spring arm.

**Note:** KL01L.IILED also can make with battery, wall type, ceiling type

### Technical Specifications

Illumination (lux) :  $\geq 60000$

Color Temperature (k) :  $4500 \pm 500$

Input power: 40 W

