

Diagnostic Equipment

Accucheck Active Glucometer

Specifications:

- Flexible dosing in or out of the meter.
- Apply blood to the test strip in the way that's easiest for you.
- User-friendly handling.
- A large display, simple handling with just 2 buttons and an easy-to-hold design.

Plus:

- Quick, 5-second results—10 seconds if dosed outside the meter
- Choice of 6 test sites on the body
- Downloadable 200-value memory with 7- and 14-day averaging
- Small 1.0 microliter sample size
- Uses 1 readily available 3-volt lithium coin cell battery
- 1,000-test battery life
- 4.6"H x 1.7"W x 0.9"D; 2.0 oz without battery
- This blood sugar test meter is compatible with ACCU-CHEK data management tools



Ecomed
Medical (Pty) Ltd

t: +27 11 955 5710
f: +27 11 955 6316

www.ecomed.co.za
sales@ecomed.co.za