

Automatic Wrist Blood Pressure Monitor

Model : S150



- Smart sense technology
- 90 memories
- One-touch automatic operation
- Carrying case



ESH approved



Smart Sense Technology



90 memories



One-touch automatic operation



AAMI SP-10

Clinically Validated

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements and this blood pressure monitor uses advanced technology that has been given the highest rating available by the prestigious Association for the Advancement of Medical Instrumentation (AAMI).

Smart Sense Technology

Smart Sense Technology utilizes an advanced air control system. The electronic deflation valve precisely controls the deflation rate and allows the monitor to detect and analyze pulse signals, adjusting cuff pressure during deflation. This significantly ensures more precise analysis, while reaching greater accuracy.



Cuff fits wrist circumference:
13.5~22cm (5.3" ~8.7")

- Storage and Transportation Condition
-10°C ~ 60°C (14°F ~ 140°F)
10% ~ 90% RH