

**ZVFF175AC VACCINE STORAGE REFRIGERATOR AND ICE PACK
FREEZER ELECTRIC – COMPRESSION**



ZERO
A P P L I A N C E S

International Cold Chain Compliant	YES
ISO 9001 : 2015	YES
IEC 60335 Compliant	YES

Specifications

Manufacturers Gross Volume (L)	Refrigerator 140L
	Freezer 32L
Hot Zone Appliance (+43°C)	YES
Freeze Protection	Grade A
Shelves	Stainless Steel
Inner Body	Stainless Steel
Fridge Temperature Display	YES
Freezer Temperature Display	YES



Performance Information

Hold Over Time During Power Cut	>48 Hours
Ice Pack Freezing	8Kg/24Hrs

Dimensions

Boxed Dimensions (mm)	(H) 1630mm x (W) 675mm x (D) 660mm
Unboxed Dimensions (mm)	(H) 1600mm x (W) 630mm x (D) 620mm
Weight (Kg's) - boxed	100kg

Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 • Int ☎ +27 10 207 1600 • Int 📠 +27 10 207 1675

🌐 info@zeroappliances.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted