

## AVC-46 (Vaccine Carrier – long range)



Model Reference	AVC-46
Type	Long Range Carrier
Hot Zone	YES
Cold Chain Compliant	YES

### Specifications

Vaccine Storage Capacity (L)	2.5
Weight - fully loaded	6.36kg
Weight - empty (with empty icepacks)	2.98kg
External Surface Material	HDPE (High Density PE)
Internal Lining Material	HIPS (High Impact PS)
Insulation Material	CFC free PU
Insulation Thickness	32/35mm

Lid Type & Fixings	Removable
Number of Ice Packs Required	4
Number of Ice Packs Supplied	4
Ice Pack Type	E5/IP.1

### Performance (E04/VC01)

Cold Life without opening	50 Hours @ 43°C
---------------------------	-----------------

### Dimensions

External Dimensions (mm)	(L) 270mm x (W) 270mm x (H) 320mm
Internal Dimensions (mm)	(L) 122mm x (W) 122mm x (H) 190mm
Vaccine Storage Dimensions (mm) with frozen ice pack	(L) 113.8mm x (W) 113.8mm x (H) 190mm



Distributed By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 • Int ☎ +27 10 207 1600 • Int 📠 +27 10 207 1675

📧 [info@zeroappliances.co.za](mailto:info@zeroappliances.co.za) • 🌐 [www.zeroappliances.co.za](http://www.zeroappliances.co.za)

E&OE Errors and Omissions Excepted