

MED85R (Single Door Medical Fridge)



Application	Fridge
Door Types	Solid Door (SD) Glass Door (GD)
ISO 9001	YES
IEC 60335 Compliant	YES

Specifications

Manufacturers Nett Volume (L)	75L
Manufacturers Gross Volume (L)	85L
Adjustable Shelves	3
Inner Body Material	Aluminium
Exterior Body Material	Stainless Steel & White Coated Mild Steel (Chromadek)
Fridge Temperature Display	YES
Temperature (Fridge)	+2°C To +8°C At +35°C Ambient
Accurate Temperature Control	YES (0.1°C Accuracy)
Visible Digital Temperature Control	YES
Internal Compartment Fan	YES
Door Locking System	YES
Glass Door Specification	Double Glazed With UV Protection (Low-E)
Automatic Defrost	YES
High / Low Temperature Alarm +2°C to +8°C	YES (Audible & Visual)
Door Open Alarm (Audible)	YES (Audible Alarm Will Sound After 60 Seconds When The Door Has Been Left Open)
Battery Alarm For Power Outage	NO (Optional)
Door Open Alarm (Visual)	YES (Visual After 60 Seconds)
CFC & HCFC Free	YES
Mains Cable Length With Plug	2.3M
Adjustable Feet	YES
Swivel Castors	NO (Optional, Add 100mm To Height)
Voltage Stabiliser	YES (Combined Surge Protect & Voltage Stabiliser)
RTMD Ready	YES
Interior Light Manual Switch	YES (External)

HOLD OVER UNIT



Performance Information

Hold Over Time During Power Cut	>48 Hours
---------------------------------	-----------

Dimensions

Internal	(H) 530mm x (W) 420mm x (D) 385mm
Boxed Dimensions	(H) 985mm x (W) 640mm x (D) 605mm
Unboxed Dimensions	(H) 975mm x (W) 620mm x (D) 585mm
Weight (Kg's) - Boxed	65



Manufactured By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

🌐 albert@zeromedical.co.za • 🌐 www.zeroappliances.co.za

E&OE Errors and Omissions Excepted