

# ECF 220 MK2 (Low Temperature Compressor Electric Freezer or Fridge)

medical refrigeration  
**MEDCO**



Max Freezer Mode Temperature	<b>-40°C</b>
Model Reference	ECF 220 MK2
Doors	1 Solid Lid

## Specifications

*Insulation Wall Thickness	120mm Polyurethane
Can Operate At +2°C To +8°C	YES
Can Operate At -40°C	YES
Gross Capacity (L)	85
Baskets Supplied	1 Top + 1 Bottom = 2
Castors	YES (4)
Lockable Lid	YES
High Temperature Alarm	YES
RTMD Ready	Optional
Voltage Protection Function	YES
ASP or AVR	Optional
Electronic Thermostat Control	YES (with display)
Climate Class (sub-tropical)	ST (up To +38°C)
Aluminium Inner Bin	YES (stainless steel optional)
Energy Sources **	Mains Electric AC 220-240V/50Hz
Colour Options	White

## Dimensions

Boxed Dimensions (mm)	(H) 970mm x (W) 825mm x (D) 680mm
Unboxed Dimensions (mm)	(H) 950mm x (W) 800mm x (D) 660mm
Weight (Kg's) - boxed	65kg

## Container Quantities

20'	42
40'	84
40'HC	84



Manufactured By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 •

📧 [albert@zeromedical.co.za](mailto:albert@zeromedical.co.za) • 🌐 [www.zeroappliances.co.za](http://www.zeroappliances.co.za)

E&OE Errors and Omissions Excepted