

## ACB-264SL (Vaccine Carrier – long range)



Model Reference	ACB-264SL
Type	Long Range Carrier
Hot Zone	YES
Cold Chain Compliant	YES

### Specifications

Gross Storage Capacity (L)	11.6
Vaccine Storage Capacity (L)	6
Weight - fully loaded	26.6kg
Weight - empty (with empty icepacks)	14.21kg
Weight Boxed - empty (with empty icepacks)	16kg
External Surface Material	LLDPE*
Internal Lining Material	LLDPE*
Insulation Material	CFC free PU
Insulation Thickness	100mm
Lid Type & Fixings	Hinged
Number of Ice Packs Required	26
Number of Ice Packs Supplied	26
Ice Pack Type	E5/IP.2
Volume of Ice Pack (L)	0.40

### Performance (E04/CB01)

Cold Life without opening	132 Hours @ 43°C
---------------------------	------------------

### Dimensions

Boxed Dimensions (mm)	(L) 650mm x (W) 530mm x (H) 460mm
External Dimensions (mm)	(L) 630mm x (W) 510mm x (H) 450mm
Internal Dimensions (mm)	(L) 400mm x (W) 280mm x (H) 206mm
Vaccine Storage Dimensions (mm)	(L) 330mm x (W) 210mm x (H) 100mm



Distributed By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📠 010 207 1675 • Int ☎ +27 10 207 1600 • Int 📠 +27 10 207 1675

🌐 info@zeroappliances.co.za • 🌐 [www.zeroappliances.co.za](http://www.zeroappliances.co.za)

E&OE Errors and Omissions Excepted