

# ECF 230 MK2 (Low Temperature Compressor Electric Freezer)

Max Freezer Mode Temperature	<b>-40°C</b>
Model Reference	ECF 230 MK2
Doors	1 Solid Lid

## Specifications

*Insulation Wall Thickness	100mm Polyurethane
Can Operate At 0°C To -40°C	YES
Gross Freezer Capacity (L)	230
Baskets Supplied	3 Top + 3 Bottom = 6
Castors	YES (4)
Lockable Lid	YES
High Temperature Alarm	YES
RTMD Ready	Optional
Voltage Protection Function	YES
ASP or AVR	Optional
Electronic Thermostat Control	YES (with display)
Temperature Zone	Up To <b>+45°C</b> (hot zone / tropical)
Aluminium Inner Bin	YES (stainless steel optional)
Energy Sources **	Mains Electric AC 220-240V/50Hz
Colour	White

## Dimensions

Boxed Dimensions (mm)	(H) 1055mm x (W) 1360mm x (D) 715mm
Unboxed Dimensions (mm)	(H) 1045mm x (W) 1330mm x (D) 710mm
Weight (Kg's) - boxed	95kg

## Container Quantities

20'	27
40'	56
40'HC	56



medical refrigeration

**MEDCO**



Manufactured By Zero Appliances (Pty) Ltd • Johannesburg, South Africa

☎ 010 207 1600 • 📞 010 207 1675 •

✉ albert@zeromedical.co.za • 🌐 [www.zeroappliances.co.za](http://www.zeroappliances.co.za)

E&OE Errors and Omissions Excepted